INSTRUCTIONS FOR USE

INTENDED USE

CD-Chex CD34® is a stabilized preparation of human blood to be used as a complete process control when evaluating CD34 positive cells. It is intended to be used with BD Biosciences ProCOUNT Progenitor Cell Enumeration Kit, BD Stem Cell Enumeration Kit, Beckman Coulter Stem-Kit, and with systems using the ISHAGE protocol for enumeration of CD34 cells.

SUMMARY AND PRINCIPLES

CD34 enumeration by flow cytometry provides a rapid and accurate assessment of the frequency of CD34 positive progenitor cells in samples from bone marrow, cord blood or peripheral blood from patients treated with hematopoietic growth factors. The ability to quantitate CD34 cells is useful in hematopoietic transplantation.

CD34 positive progenitor cells can be distinguished on the basis of light scatter properties in conjunction with surface antigens. CD-Chex CD34 is designed to represent a blood sample containing CD34 positive cells having characteristics similar to progenitor cells: low/intermediate side scatter properties, CD34 expression and low expression of CD45 (compared to lymphocytes).^{1,2,3} When stained with monoclonal antibodies for CD34 positive cell enumeration, CD-Chex CD34 control will provide reference values for CD34 positive cells within the ranges on the assay sheet.

REAGENTS

CD-Chex CD34 contains human leukocytes and erythrocytes in a preservative medium.

PRECAUTIONS

- For In Vitro Diagnostic Use.
- CAUTION:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at streck.com for specific FDA required blood tests.
- This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
- This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates any diagnostic use of the product.
- CD-Chex CD34 products should not be used as a calibrator.
- SDS can be obtained at streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

CD-Chex CD34 is stable through the expiration date when stored at 2 °C to 10 °C. After opening, CD-Chex CD34 is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C. DO NOT FREEZE.

INDICATIONS OF PRODUCT DETERIORATION

If CD-Chex CD34 values are not within the expected range on the assay sheet or granulocytes show a loss of forward scatter (FSC):

- Review control product package insert and the operating procedure of the instrument.
- Check expiration date of the product on the vial. Discard outdated products.
- Assay an unopened vial of CD-Chex CD34. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.
- Clumping of the cell suspension indicates instability or deterioration, in which case the product should not be used.

INSTRUCTIONS FOR USE

CD-Chex CD34 is designed to be used with ISHAGE CD34 enumeration protocols.^{1,2,3}

- Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
- Remove a vial of the control from refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
- Mixing Procedure (**mechanical mixing by vortex or rotator is not recommended**):

For a video demonstration, visit streck.com/mixing.

- Holding the vial vertically between the palms of the hands, roll the vial back and forth for 20-30 seconds.

- Hold the vial by the ends between the thumb and finger, and mix by gently inverting the vial at least 8-10 times end-over-end until all cells are thoroughly suspended.

- Aliquot immediately after mixing.
- Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to sampling.

Note: Vials stored for an extended period of time may require extra mixing.

- Return control reagent to refrigeration after sampling to ensure maximum open-vial stability.
- Add recommended monoclonal antibodies according to manufacturer's instructions to each tube and mix gently.
- Incubate according to antibody manufacturer's instructions.
- Add recommended amount of RBC lysing agent and follow manufacturer's instructions.
- Analyze by flow cytometry using your laboratory's established protocol.

LIMITATIONS

- CD-Chex CD34 is not recommended for use with antibodies that target the class I and class II epitopes of CD34.
- CD-Chex CD34 is designed to be used with a RBC lysing agent and may not provide results within the assay range if analyzed without RBC lysis.

EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on properly compensated flow cytometers.^{1,2,3,4} The assay values are obtained using common flow cytometry reagents. See the IFU and assay for limitations or specific instructions for reagents.

Upon receipt of a new control lot, it is recommended that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different reagents, laboratory protocols, instrument calibration, maintenance, and operator technique.

REFERENCES

- Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
- Keeney M, Chin-Yee I, Weir K, Popma J, Nayar R, Sutherland DR. Single platform flow cytometric absolute CD34+ cell counts based on ISHAGE guidelines. *Cytometry (Communications in Clinical Cytometry)* 1998; 34:61-70.
- Sutherland DR, Anderson L, Keeney M, Nayar R, Chin-Yee IH. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. *Journal of Hematotherapy* 1996; 5:213-226.
- Marti G, Johnsen H, Sutherland R, Serke S. Letter to the editor: A convergence of methods for a worldwide standard for CD34+ cell enumeration. *Journal of Hematotherapy* 1998; 7:105-109.

QUALITY CONTROL PROGRAM

Streck offers *STATS*®, an interlaboratory quality control program, to all customers at no charge. For more information, contact the *STATS* Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at streck.com.

See streck.com/patents for patents that may be applicable to this product.

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LITTERATUR

1. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
2. Keeney M, Chin-Yee I, Weir K, Popma J, Nayar R, Sutherland DR. Single platform flow cytometric absolute CD34+ cell counts based on ISHAGE guidelines. Cytometry (Communications in Clinical Cytometry) 1998; 34:61-70.
3. Sutherland DR, Anderson L, Keeney M, Nayar R, Chin-Yee IH. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. Journal of Hematotherapy 1996; 5:213-226.
4. Marti G, Johnsen H, Sutherland R, Serke S. Letter to the editor: A convergence of methods for a worldwide standard for CD34+ cell enumeration. Journal of Hematotherapy 1998; 7:105-109.

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RÉFÉRENCES

1. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
2. Keeney M, Chin-Yee I, Weir K, Popma J, Nayar R, Sutherland DR. Single platform flow cytometric absolute CD34+ cell counts based on ISHAGE guidelines. *Cytometry (Communications in Clinical Cytometry)* 1998; 34:61-70.
3. Sutherland DR, Anderson L, Keeney M, Nayar R, Chin-Yee IH. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. *Journal of Hematotherapy* 1996; 5:213-226.
4. Marti G, Johnsen H, Sutherland R, Serke S. Letter to the editor: A convergence of methods for a worldwide standard for CD34+ cell enumeration. *Journal of Hematotherapy* 1998; 7:105-109.

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BIBLIOGRAFIA

- Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.
- Keeney M, Chin-Yee I, Weir K, Popma J, Nayar R, Sutherland DR. Single platform flow cytometric absolute CD34+ cell counts based on ISHAGE guidelines. Cytometry (Communications in Clinical Cytometry) 1998; 34:61-70.
- Sutherland DR, Anderson L, Keeney M, Nayar R, Chin-Yee IH. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. Journal of Hematotherapy 1996; 5:213-226.
- Marti G, Johnsen H, Sutherland R, Serke S. Letter to the editor: A convergence of methods for a worldwide standard for CD34+ cell enumeration. Journal of Hematotherapy 1998; 7:105-109.

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