## INSTRUCTIONS FOR USE INTENDED USE

Para 12® Extend is an assayed hematology control for evaluating the accuracy and precision of hematology instruments that provide a white blood cell differential.

## SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated, semi-automated and manual procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation.

Use of stabilized cell preparations for controlling hematology instrumentation is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument, Para 12 Extend will provide values within the expected range indicated on the assay sheet.

## REAGENTS

This whole blood reagent may contain any or all of the following: stabilized human red blood cells, a white blood cell component consisting of human and/or non-human analogs and a platelet component consisting of a human or non-human analog in a preservative medium.

## PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at [www.streck.com](http://www.streck.com) for specific FDA required blood tests.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
5. Control products are not to be used as calibrators.
6. SDS can be obtained at [www.streck.com](http://www.streck.com), by calling 800-843-0912, or by calling your local supplier.

## STORAGE AND STABILITY

Para 12 Extend is stable through the expiration date when stored at 2 °C to 10 °C. After opening, Para 12 Extend is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C.

## INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of Para 12 Extend. Discard outdated products.
3. Assay an unopened vial of Para 12 Extend. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or [technicalservices@streck.com](mailto:technicalservices@streck.com).


## INSTRUCTIONS FOR ASSAY ZIP FILE DOWNLOAD

### Windows® operating systems

1. Insert a formatted USB flash drive or a formatted 3½-inch floppy disk into the computer.
2. Go to the product page at [www.streck.com](http://www.streck.com) and download the assay zip file for the lot number you are using.
3. When prompted, click "Open" to download and open the file.
4. Once the file is opened, click "File" → "Extract All" or "Extract All..."  
Note: Windows 7 users - Be sure the box reading "Show extracted files when complete" is unchecked.
5. Follow the "Extraction Wizard" (if applicable). Click on the "Browse" button to select an extract destination. Select the flash drive or floppy disk that was inserted in step 1. Click extract.
6. After the files have been extracted, safely remove the flash drive or floppy disk from the computer and insert into the appropriate instrument.
7. Refer to instrument operator's manual for assay file set up.

## INSTRUCTIONS FOR USE

1. Remove vials of control from the refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
2. To mix: **(Do NOT mix mechanically or vortex.)**  
**For a video demonstration, visit [www.streck.com/mixing](http://www.streck.com/mixing).**
  - a. Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- b. Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- c. Analyze immediately after mixing. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to instrument analysis.

- d. Steps a-c must be repeated upon removing the sample from the refrigerator for the entire open-vial time period regardless of the method of analysis (open tube, cap piercing, auto sample or manual sample).
3. Refer to the instrument manual for the system in use for analyzing control materials.
  4. After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both the vial and cap before replacing cap and returning to refrigeration.

## LIMITATIONS

A manual differential analysis of white blood cells cannot be accomplished with Para 12 Extend. The white blood cell components simulate white blood cells in size, not morphology.

## EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments as well as by manual reference methods. MCV and hematocrit values are not corrected for trapped plasma. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique.

Assay values and expected ranges for instruments not listed on the assay sheet must be established by the user. Per CLSI H26-A2<sup>1</sup> it is recommended that each level of control be run twice a day for 3-5 days to establish individual lab means for each measurand.

Trending in the MCV parameter over the product shelf life is inherent to hematology control products. This characteristic does not indicate product instability. Ranges and limits for these parameters may need to be adjusted.

## REFERENCES

1. Clinical and Laboratory Standards Institute, H26-A2, Validation, verification, and quality assurance of automated hematology analyzers; Approved Standard - Second Edition.

## QUALITY CONTROL PROGRAM

Streck offers *STATS*®, an interlaboratory quality control program, to all customers at no charge. For more information, contact the *STATS* Department at 800-898-9563 or [statsdata@streck.com](mailto:statsdata@streck.com). Additional information can be found at [www.streck.com](http://www.streck.com).

## ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at [www.streck.com](http://www.streck.com).

## Rx Only

## GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at [www.streck.com](http://www.streck.com).

See [www.streck.com/patents](http://www.streck.com/patents) for patents that may be applicable to this product.

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