INSTRUCTIONS FOR USE

INTENDED USE

STaK-Chex[®] Plus Retics is an assayed whole blood control for evaluating the accuracy and precision of automated procedures that measure blood cell parameters.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated procedures that measure whole blood parameters. Daily use of this whole blood control provides quality control data for confirming the precision and accuracy of instrument operation protocols.

Use of stabilized cell preparations for controlling laboratory testing protocols is an established procedure. When handled like a patient sample and assayed on a properly calibrated and functioning instrument or manual methods, this whole blood control will provide values within the expected range indicated on the assay sheet.

REAGENTS

These whole blood reagents may contain any or all of the following: stabilized human red blood cells, a nucleated red blood cell analog, a white blood cell component consisting of human or non-human analogs and a platelet component consisting of a human or non-human analog in a preservative medium.

PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at www.streck.com for specific FDA required blood tests.
3. This product should not be disposed of in general waste, but should be disposed of with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates any diagnostic use of the product.
5. Control products are not to be used as calibrators.
6. SDS can be obtained at www.streck.com, by calling 800-843-0912, or by calling your local supplier.

STORAGE AND STABILITY

STaK-Chex Plus Retics is stable through the expiration date when stored at 2 °C-10 °C before and after opening. After opening, STaK-Chex Plus Retics is stable throughout the open-vial expiration dating, as indicated on the assay sheet, when stored at 2 °C-10 °C.

INDICATION OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. Discoloration of the product may be caused by overheating or freezing during shipping or storage. Darkly colored supernatant may be indicative of product deterioration, however, moderately colored supernatant is normal and should not be confused with product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of the product on the vial. Discard outdated products.
3. Assay an unopened vial of the control. If the values are still outside the Expected Range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.

INSTRUCTIONS FOR USE

1. Remove vials of STaK-Chex Plus Retics from the refrigerator and warm to room temperature (18 °C-30 °C) for 15 minutes before use.
2. To mix: **(Do NOT mix mechanically or vortex.)**
For a video demonstration, visit www.streck.com/mixing.
a. Hold the vial vertically and roll each vial between the palms of the hands for 15-20 seconds.

- b. Continue to mix by holding the vial by the ends between the thumb and finger, rapidly inverting the vial 20 times end-over-end using a very quick turning motion of the wrist.

- c. Analyze immediately after mixing. Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to instrument analysis.
 - d. Steps a-c must be repeated upon removing the sample from the refrigerator for the entire open-vial time period regardless of the method of analysis (open tube, cap piercing, auto sample or manual sample).
3. Refer to the instrument manual for analyzing control material.
 4. After sampling, return to refrigeration for maximum open-vial stability. If run in the open mode, wipe the threads of both the vial and cap before replacing cap and returning to refrigeration.

LIMITATIONS

A manual differential analysis of white blood cells cannot be accomplished with this product. The white blood cell components simulate white blood cells in size, not morphology.

EXPECTED RESULTS

The mean assay values provided for each parameter are derived from replicate analyses on calibrated instruments as well as by manual reference methods. If the product is value assigned for MCV and hematocrit, assay values are not corrected for trapped plasma. The assay values are obtained using reagents recommended by instrument manufacturers and are to be used for instrument control; they are not absolute assays for calibration. Upon receipt of a new control lot, it is good laboratory practice

that an individual laboratory establish its own mean and limits for each parameter. However, the control means established by the laboratory should fall within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance, and operator technique.

Assay values and expected ranges for instruments not listed on the assay sheet must be established by the user. Per CLSI H26-A2, it is recommended that each level of control be run twice a day for 3-5 days to establish individual lab means for each measurand.

'R' alarms may occur with control materials on automated systems. This will not affect the validity of the results. The alarms may be disregarded if the control is performing within the assay ranges.

Trending in the MCV parameter over the product shelf life is inherent to hematology control products. This characteristic does not indicate product instability. Ranges and limits for these parameters may need to be adjusted.

ONLINE ASSAY FILE DOWNLOAD INSTRUCTIONS

For Beckman Coulter[®] LH780/785, LH750/755, LH500 & HmX

1. Insert a blank USB flash drive or 3½-inch floppy disk into the computer.
2. On the appropriate product page at streck.com, click the link for the assay disk file with the corresponding lot number you are using.
3. When prompted, Click "Open" to download and open the file.
4. Once the file is opened, click "File" → "Extract All" or "Extract All..."
Note: Windows[®] 7 users - Be sure the box reading "Show extracted files when complete" is unchecked.
5. Follow the "Extraction Wizard" (if applicable). Click on the "Browse" button to select an extract destination. Select the flash drive or floppy disk that was inserted in Step 1.
6. After the files have been extracted, safely remove the flash drive or floppy disk from the computer and insert into the appropriate instrument.
7. Continue following appropriate control file set up procedures.

REFERENCES

1. Clinical Laboratory Standards Institute, H26-A2, Validation, verification, and quality assurance of automated hematology analyzers. Approved Standard - Second Edition.

QUALITY CONTROL PROGRAM

Streck offers STaTS[®], an interlaboratory quality control program, to all customers at no charge. For more information, contact the STaTS Department at 800-898-9563 or statsdata@streck.com. Additional Information can be found at www.streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

Rx Only

GLOSSARY OF SYMBOLS

See the Instructions (IFU) tab under Resources on the product page at www.streck.com.

The brand and product names of the instruments are trademarks of their respective holders.

See www.streck.com/patents for patents that may be applicable to this product.

Streck
7002 S. 109 Street, La Vista, NE 68128 USA

MEDI MARK[®] Europe
11, rue Emile Zola, BP 2332
38033 Grenoble Cedex 2, France

350289-22
2019-09

CONTROL FILE SETUP PROCEDURES FOR BECKMAN COULTER® LH500, LH750/755, LH780/785; HmX; UniCel® DxH™ 600 & UniCel DxH 800

LH 500 - QC SETUP

Separate CBC/Diff and Retic files

1. Insert the USB flash drive or 3½-inch floppy disk with the assay values for the new lot.
 2. Select the yellow diamond icon \diamond on the Command Center to display the System Setup window.
 3. Select the QA icon, select the Controls tab and click on the icon, New Control folder \mathbb{E} .
 4. Select the following:
 - a. Source - **BCI/Coulter**
 - b. Type - **5C Series**
 - c. Level - **All levels**
 5. Select the icon Set Up New Lot to display the Diskette Entry screen.
 6. Select Auto Transmit and click on the OK icon. The assay information will be transmitted to the QC CBC/Diff files.
 7. To download the Retic information, repeat the steps above, entering **Retic C** in the Type field (step 4b).
- Low values will load into the CBC Abnormal I file. (Retic Level I)**
Normal values will load into the CBC Normal file. (Retic Level II)
High values will load into the CBC Abnormal II file. (Retic Level III)

LH 500 - RUNNING CONTROLS

CBC/Diff file

1. In the Command Center located at the bottom of the screen, select the following:
 - a. Process Type - **Auto Analysis**
 - b. Run Type - **CD (CBC/Diff)**
 - c. Asp Mode - **Auto**
2. Position the vials in the cassette holder so the instrument can scan the **CBC/Diff barcode**.
3. Place the cassette into the loading bay and the instrument will begin to process the controls.

Retic file

1. For Retic analysis, consult with instrument Operator Manual for instructions on how to prepare Retic samples and repeat the steps above entering **R** for Run Type (step 1b).
2. Enter the Retic lot number into the barcode field and press enter.

LH 750/755 & 780/785 - QC SETUP

Separate CBC/Diff and Retic files

1. Insert the USB flash drive or 3½-inch floppy disk with the assay values for the new lot.
 2. Select the yellow diamond icon \diamond on the Command Center to display the System Setup window.
 3. Select the QA icon, select the Controls tab and click on the New Control folder icon \mathbb{E} .
 4. Select the following:
 - a. Source - **BCI/Coulter**
 - b. Type - **5C**
 - c. Level - **All levels**
 5. Select the Set Up New Lot icon to display the Diskette Entry screen.
 6. Select Auto Transmit and click on the OK icon. The assay information will be transmitted to the QC CBC/Diff files.
 7. To download the Retic information, repeat the steps above, entering **Retic C** in the Type field (step 4b).
- Low values will load into the CBC Abnormal I file. (Retic Level I)**
Normal values will load into the CBC Normal file. (Retic Level II)
High values will load into the CBC Abnormal II file. (Retic Level III)

LH 750/755 & LH 780/785 - RUNNING CONTROLS

Separate CBC/Diff and Retic files

1. NRBC parameters can be recovered by enabling auto-print for the controls and manually recording the values.
2. In the Command Center located at the bottom of the screen, select the following:
 - a. Process Type - **Auto Analysis**
 - b. Default Type - **CD (CBC/Diff)**
 - c. Random Access - **Select the box**
3. Position the vials in the cassette holder so the instrument can scan the **CBC/Diff barcode**.
4. Place the cassette into the loading bay and the instrument will begin to process the controls.
5. When the CBC/Diff analysis is complete, repeat the steps above for the Retic analysis entering **R** for Default Type (step 2b):
6. Position the vials in the cassette holder so the instrument can scan the **Retic barcode**.
7. Place the cassette into the loading bay and the instrument will begin to process the controls.

LH 750/755 & LH 780/785 - QC SETUP

Combined CBC/Diff/Retic file

1. Select the yellow diamond icon \diamond on the Command Center to display the System Setup window.
2. Select the QA icon, select the Controls tab, and click on the icon for New Control folder \mathbb{E} .
3. Select the following:
 - a. Source - **Other**
 - b. Type - **CBC/Diff/Retic**
 - c. Level - **Low, Normal or High**
4. Select the icon, Set Up New Lot, to display the control setup screen.
5. Enter the 8-digit CBC/Diff/Retic lot number found in the lot number chart on the assay. (ex. 03314187 - Low)
6. Enter the expiration date, assay means and expected ranges from the assay for each level. Select the OK icon to save the data.
7. Repeat steps 3 - 5 for the remaining levels. All values will be transmitted into one control file for each level.

LH 750/755 & LH 780/785 - RUNNING CONTROLS

Combined CBC/Diff/Retic file

1. In the Command Center located at the bottom of the screen, select the following:
 - a. Process Type - **Auto Analysis**
 - b. Default Type - **CDR (CBC/Diff/Retic)**
 - c. Random Access - **Select the box**
2. Position the vials in the cassette holder so the instrument can scan the **CBC/Diff/Retic barcode**.
3. Place the cassette into the loading bay and the instrument will begin to process the controls.
4. **Note: NRBC parameters can be recovered by enabling auto-print for the controls and manually recording the values.**
 - a. Select Run Configuration in the Command center, and select Control in the QA Samples.
 - b. Run controls and record NRBC values from printout.

HmX - QC SETUP

Expected Range file

1. From the Control Setup screen, select Expected Range and the appropriate file type and level.
2. Enter the Expected Ranges for each parameter as listed on the assay sheet. For CBC/Diff files:
Enter the Low (87) ranges into the Abnormal I file.
Enter the Normal (88) ranges into the Normal file.
Enter the High (86) ranges into the Abnormal II file.

For Retic files:

- Enter the Low (31) ranges into the Retic Level I file.
- Enter the Normal (42) ranges into the Retic Level II file.
- Enter the High (53) ranges into the Retic Level III file.

3. Press F10/ESC to save.

Note: The Expected Ranges only need to be entered once as they do not change from lot to lot.

HmX - See your instrument operator's manual for instructions on how to upload assay values. HmX retic assay values must be entered manually.

CBC/Diff file

1. From the Control Setup screen, select **CBC/Diff Files**.
2. Enter the 6-digit lot number and expiration date from the assay sheet. The expected ranges for each level will be transmitted into QC files after the lot number is entered.
3. Enter other applicable information for your laboratory.
4. Enter the assigned values from the assay sheet.
5. Press F10 Save/ESC to save.

Retic file

1. From the Control Setup screen, select **Retic Files**.
2. Enter the 6-digit Retic lot number located near the bottom of the assay sheet. The expected ranges for each level will be transmitted into the QC files once the lot number is entered.
3. Enter the expiration date and other applicable information for your laboratory.
4. Enter the Ret% value from the assay sheet. In order to have the Ret# calculated, enter the RBC value from the assay in the Ref RBC field.
5. Press F10 Save/ESC to save.

HmX - RUNNING CONTROLS

CBC/Diff file

1. Position the vials in the cassette holder so the instrument can scan the **CBC/Diff barcode**.
2. Place the cassette into the loading bay and the instrument will begin to cycle the controls.

Retic file

1. Prepare a vial and change the instrument mode to analyze reticulocytes.
2. On the numeric keypad, press ID and enter the 6-digit Retic lot number.
3. Aspirate the prepared sample.

UniCel DxH 600/800 - QC SETUP (Presence of RUO parameters designated with the @ symbol will not allow assay upload using the barcodes.)

Combined CBC/Diff/Retic files

1. Select: **Menu/Setup/Quality Control**
2. Select: **New Patient Control**
3. Confirm the cursor selection is in the Specimen ID field: **Scan lot number barcode**
4. Choose Type: **CDR**
5. Choose Level
6. Select: **Create Control**
7. Confirm the cursor selection is in the Control ID field: **Scan 2-D barcode**
8. Select: **Save**
9. Repeat steps 2-8 for additional Control Levels

UniCel DxH 600/800 - RUNNING CONTROLS

Combined CBC/Diff/Retic

1. Place STKR vials in cassette with CBC/Diff/Retic barcode facing out.
2. Place cassette into output buffer
3. Data will automatically transfer to appropriate "Patient Control" file.

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