

### Instruction for Use

# **NukEx PLUS**

## **Nucleic Acid Release Reagent**

Reagent for the enzymatic release of DNA from bacterial cultures and swabs.



#### G01073



100





gerbion GmbH & Co. KG

Remsstr. 1

70806 Kornwestheim

Germany

phone: +49 7154 806 20 0 fax: +49 7154 806 20 29 e-mail: info@gerbion.com

www.gerbion.com

## Index

1	Intended Use	. 4
2	Principle of the Test	. 4
3	Package Contents	. 4
4	Equipment and Reagents to be supplied by User	. 4
5	Transport and Storage	. 4
6	Important Notes	. 4
7	Waste Handling	. 5
8	Sample Material	. 5
9	Procedure	. 5
10	Storage of Crude Lysates	. 5
11	Troubleshooting	. 6
12	Abbreviations and Symbols	. 6
13	Literature	F

#### 1 Intended Use

The NukEx PLUS Nucleic Acid Release Reagent is designed for the release of DNA from bacterial cultures and swabs. The crude lysates can be directly applied in molecular assays, e.g. real time PCR.

#### 2 Principle of the Test

The NukEx PLUS Nucleic Acid Release Reagent causes the lysis of bacteria. The lysis is performed for 15 minutes at room temperature.

Nucleic acids released with NukEx PLUS Nucleic Acid Release Reagent can be analysed by employing the supernatants obtained directly in the subsequent molecular assay.

Pooling of the lysates prior to analysis is possible; however, it is subject to the purpose and regulations of the particular application.

#### 3 Package Contents

4x 8.75 ml NukEx PLUS Nucleic Acid Release Reagent sufficient for 100 reactions.

#### 4 Equipment and Reagents to be supplied by User

- Laboratory equipment according to national safety instructions
- Sterile pipet tips with filter
- Nuclease-free 1.5 or 2.0 ml microcentrifuge tube
- Optional: Liquid handling systems for automation (e.g. Myra, Bio Molecular Systems)

#### 5 Transport and Storage

The NukEx PLUS Nucleic Acid Release Reagent must be stored  $\leq$  -18°C. If properly stored, the product is stable until the date of expiry printed on the label.

Please note, that improper storage will adversely impact nucleic acid release due to decreased enzymatic activity. Therefore, NukEx PLUS Nucleic Acid Release Reagent is always shipped on dry ice.

#### 6 Important Notes

- The Nukex PLUS Nucleic Acid Release Reagent must be utilized by qualified personnel only.
- Good Laboratory Practice (GLP) has to be applied.

 Clinical samples must always be regarded as potentially infectious material and all equipment used has to be treated as potentially contaminated.

#### 7 Waste Handling

- Dispose of unused reagents and waste should occur in accordance with country, federal state and local regulations.
- Material Safety Data Sheets (MSDS) are available upon request.

#### 8 Sample Material

Starting material are respiratory swabs (e.g. buccal swabs, nasal swabs, etc.) or bacterial culture samples.

#### 9 Procedure

- Pipet 350 μl NukEx PLUS Nucleic Acid Release Reagent into an appropriate tube (e.g. 2 ml reaction tube, safe lock).
- Place the swab tip or the picked bacterial colony into the reaction tube and break or cut off the applicator at a length that allows the tube to be closed.
- Close reaction tube tightly.
- Vortex thoroughly 4- 5 times
- Incubate for 15 min at room temperature.

#### 10 Storage of Crude Lysates

For storage conditions of crude NukEx PLUS Nucleic Acid Release Reagent lysates please refer to table 1.

Table 1: Storage conditions for crude lysates

Time	Storage Condition	
up to 6 hours	Room temperature	
up to 24 hours	+2 to +8 °C	
long term storage	≤ - 18°C	

#### 11 Troubleshooting

The following troubleshooting guide is included to help you with possible problems that may arise in a subsequent PCR.

Neither sample nor Internal Control show a PCR signal						
Concentration of PCR	Components present in the sample may inhibit					
inhibitors in the sample	the PCR. Therefore, dilute the supernatant 1:10					
too high	in dH <sub>2</sub> O (PCR grade). If necessary, extract the					
	nucleic acid from the crude lysate with a					
	commercial extraction kit (e.g. NukEx Pure					
	RNA/DNA) and repeat PCR analysis.					
Taq-Polymerase	PCR run should be started immediately (max.					
damaged by NukEx	within 15 – 20 min.) after pipetting the NukEx					
PLUS	PLUS lysate into the reaction mix.					
Negative PCR result for a known-positive sample, Internal Control shows						
no inhibition						
Kit stored under non-	Store kit at ≤ -18°C. Do not use after the date of					
optimal conditions or	expiry printed on the label.					
kit expired						
Incorrect incubation	Make sure incubation conditions comply with					
conditions	the protocol.					

#### 12 Abbreviations and Symbols

DNA	Desoxyribonucleic Acid	\$	<pre>contains sufficient for <n> test</n></pre>
PCR	Polymerase Chain Reaction	18°C	Upper limit of temperature
LOT	Batch code	***	Manufacturer
CONT	Content	$\subseteq$	Use by YYYY-MM
REF	Catalog number	i	Consult instructions for use

#### 13 <u>Literature</u>

[1] Sambrook, J. and Russell, D.W.: Molecular Cloning, 2001.