



Retic-Chex[®] Linearity



Retic-Chex Linearity is designed to establish the reportable range and linearity of select hematology analyzers that report the reticulocyte parameter. The linearity control kit is assayed for Sysmex and Beckman Coulter hematology analyzers.

Features

- + Plastic vials with pierceable caps for analyzer autosampling
- + 105-day closed-vial stability
- + Assayed for Sysmex and Beckman Coulter hematology analyzers
- + 5-day open-vial stability

Benefits

- + Helps clinical laboratories meet quality requirements
- + Established assay ranges provide less room for error
- + Access to *STATS*[®], a free interlaboratory quality control program with peer group data comparison and online submission

Ordering information

Description	Catalog#
5 x 2.5 mL for Sysmex	285532
5 x 3.0 mL for Beckman Coulter	285535